

Author(s): JA, BH

09/03/2013

Protocol for Long-Term Islet Culturing

1. PURPOSE:

This protocol describes how to culture islets for long term purposes.

2. MATERIALS REQUIRED:

1. PIM(S)[®] media (Prodo labs – Cat no. - PIM-S001GMP).
2. PIM(R)[®] media (Prodo labs-Cat no. – PIM-R001GMP).
3. PIM(ABS)[®] (Prodo labs – Cat no. - PIM-ABS001GMP).
4. PIM(G)[®] (Prodo labs – Cat no. - PIM-G001GMP).
5. 35x10mm petri dish (Fisher, Sterile (Cat no. – 08-757-100A)).
6. T-150 non-tissue culture treated flasks (Corning* Non-Treated Culture Flasks, Polystyrene, Sterile, (Cat no. - 431465).

3. PROCEDURE:

To begin, follow the **Recovering Islets After Shipping Protocol**. To avoid or minimize the chance of contamination, the appropriate steps below are to be performed in a laminar flow hood with good sterile technique.

• ISLET CULTURE

Short term islet culture is done in a 37°C incubator, with 5% CO₂, following the Culturing Islets Protocol.

1. The PIM (S)[®] media on the islets in short term culture needs to be changed every 3-4 days. This is done by following the **Protocol for Islet Media Change**.
2. Also, each time the media is changed, take islet samples for assessing viability/purity and place them in a 35mm petri dish.
3. Islets can be cultured up to 10 days using this method.

• EXTENDED CULTURING OF ISLETS (AFTER 10 DAYS)

Extended culturing is done using PIM(R) in a 37°C incubator, with 5% CO₂

1. On the 10th day following the isolation, the islets need a 50% media change to PIM(R)[®]. This is done following the **Protocol for Islet Media Change**.
2. The flasks are then placed in the 37°C incubator, with 5% CO₂.
3. After this, step 1 is repeated every 3-4 days, using PIM(R)[®]
4. Islets can be cultured for at least 4 weeks using this method